



Always Safe

Use this checklist when placing a child into equipment and you will always be safe in the knowledge that the child is secure and comfortable.

Check...



the lap belt and supporting straps are tight enough

The most important adjustment for a child with postural problems. Allow a finger's width between belt and child, depending on clothing.



any excess length on supporting straps is securely fastened

Preventing siblings or class mates undoing straps.



seat depth

It's important that the legs are fully supported, check there's a two finger width behind the child's knees to prevent chafing.



all adjustments are working and set correctly

Lock these in position to prevent siblings or class mates tampering with them!



castors are locked

When the product is stationary.



equipment height is correct

This gives your child the chance to interact with their surroundings and socialise with their friends.



all surfaces are fully clean

Wipe down if necessary.



instructions for use & service record

Always read the equipment manufacturer's instructions for use and check the service record.

Double check...

...again, all lap belts and supporting straps are securely fastened to ensure each child receives the full postural support they need.