



SleepSafe Beds, LLC  
3629 Reed Creek Drive  
Bassett, VA 24055

Local: 276-627-0088  
Fax: 276-627-0234  
SleepSafeBed.com

Dear DME Professional,

We have evaluated the compatibility of the Jenx® Dreama: 24hr Positioning System for use in SleepSafe® Beds. There are two accommodations required that ensure optimized safety and effectiveness.

- 1) One (1) Jenx® Dreama system link must be removed to adapt properly to the length of the SleepSafe Bed deck.
- 2) A specially fitted pad and strap, provided by SleepSafe® Beds, secures the Dreama: 24hr Positioning System to the bed deck. This pad provides additional stability to the system while also shielding the user from potential pinch points.

When these requirements are met, the SleepSafe® bed becomes a suitable foundation for users who require the Jenx® Dreama© 24hr Positioning System for their therapy and well-being.

Best regards,

Bob Baldacci

Director of Quality and Regulatory Affairs, *SleepSafe® Beds*

866-852-2337 | [www.SleepSafeBed.com](http://www.SleepSafeBed.com)