

BILLING/SHIPPING INFORMATION					
Date:			Name:		
PO Number:			Address:		
Branch:					
Mark for:			City:		
Dealer Contact:			State:		
Email Address:			Zip Code:		
Phone:			Phone:		
<b>QUOTE</b>			<b>ORDER</b>		

A versatile standing system offering prone or upright standing that accommodates abducted or neutral leg positioning for children between approx 1-9 years.



### KEY FEATURES

- Tool free adjustments
- Individual leg abduction and adjustment
- Anti-microbial covers to aid infection control in multi-use environments
- Built in gauges enable precise positioning

	<b>Size 1 Prone</b>
Age Range (approx)	6 months* - 9 years
User Height	21"* - 54.3"
User Weight (max)	99 Lbs
Chest/Hip Pad Width (Min/Max)	7" - 12"
Height from footplates to top of upper support	23" - 9"
Height of Stander (floor to top of upper pad, <b>upright</b> position)	43"
Height of Stander (floor to support board, <b>horizontal</b> position)	32"
Kneeblock Width (center to center, in neutral standing position)	5" - 8"
Footplate Angle Range	10 - 25 degrees
Individual Leg Abduction	Neutral to 30 degrees
Angle Range	0 - 90 degrees

\*Measurements can be achieved by the use of additional mini kit accessories (ST022)

### Suggested HCPCS code E0638\*

*\*It is the responsibility of the CRT provider to determine the HCPCS code most appropriate for each circumstance. The HCPCS codes for reference only and are not to be deemed as legal advice. Actual coding and coverage varies by funding entities and the use of HCPCS does not guarantee coverage or payment.*

**CHOOSE YOUR COLOR- All part numbers containing \*\* will color match**



## REQUIRED COMPONENTS

<b>BASE</b>		<b>Compact base</b> with gas strut assisted angle adjustment and incorporated angle gauge and safety lock. A: 35" B: 25"	<b>MAST</b>		<b>Mast with dual leg posts and individual footplates</b>	<b>SUPPORT BOARDS</b>		<b>Support Boards.</b> Pair in color of choice. Top board height adjustable to accommodate growth.
	STZ-BS1-PRN	\$1582		<input checked="" type="checkbox"/>	STZ-M1		\$985	<input checked="" type="checkbox"/>

<b>TRUNK SUPPORT</b>		<b>Lateral Supports with Rotation Control Strap</b> Width and height adjustable - can be set to accommodate asymmetry. With incorporated rotation control strap for simple and effective control and stability.
	STZ012-**	\$466 <input type="checkbox"/>

OR

	<b>Rotation Control Strap</b> A wide support belt without lateral pads, threaded through the support fixings. Fastened securely around child and then pulled to facilitate extension or rotation.
STZ013-**	\$466 <input type="checkbox"/>

<b>PELVIC SUPPORT</b>		<b>Lateral Supports with Rotation Control Strap</b> Width adjustable with incorporated rotation control strap for simple and effective control and stability.
	STZ012-**	\$466 <input type="checkbox"/>

OR

	<b>Rotation Control Strap</b> A wide support belt without lateral pads, threaded through the support fixings. Fastened securely around child and then pulled to facilitate extension or rotation.
STZ013-**	\$466 <input type="checkbox"/>

## REQUIRED COMPONENTS: continued

<b>KNEE CONTROL</b>		<p><b>Standard Knee Support</b> Width adjustable side walls mounted on an angle and depth adjustable bracket. Tool free adjustment. <b>PAIR</b></p>	<p>A: 5" B: 5" C: 3" - 4"</p>	<p><b>STZ 007-**</b> \$758</p>	<input type="checkbox"/>
		<p><b>Standard Knee Strap</b> holds knee in desired position and prevents flexion. <b>PAIR *Optional</b></p>	<p>A: 7" B: 17"</p>	<p><b>MU54-**</b> Size 1 \$134</p>	<input type="checkbox"/>
		<p><b>De-Rotational Knee Straps</b> Allows for cuffing of the child's leg to control internal or external knee rotation. <b>PAIR *Optional</b></p>	<p>A: 15" B: 12" C: 12"</p>	<p><b>MU45-**</b> Size 1 \$152</p>	<input type="checkbox"/>
		<p>A: 10" B: 21"</p>	<p><b>MU68-**</b> Size 2 \$175</p>	<input type="checkbox"/>	
		<p>A: 18" B: 11." C: 11"</p>	<p><b>MU69-**</b> Size 2 \$175</p>	<input type="checkbox"/>	

<b>FOOT POSITIONING</b>		<p><b>Sandals</b> Ensure solid foot placement for effective leg positioning and standing. <b>Pair</b></p>	<p>Size 1- A: 3" B: 7"</p>	<p><b>S16</b> \$178</p>	<input type="checkbox"/>
		<p>Size 2- A: 4" B: 9"</p>	<p><b>S2</b> \$205</p>	<input type="checkbox"/>	

OR

	<p><b>Multi-Axis Knee Supports-</b> Includes Arched Foam Pad <b>and</b> support straps mounted on multi-axis bracket. <b>PAIR</b></p>	<p>A: 5" B: 5" Max knee circumference: 12"</p>	<p><b>STZ017-**</b> \$1072</p>	<input type="checkbox"/>
	<p>A: 5.3" B: 4.5" Max knee circumference: 18"</p>	<p><b>STZ023-**</b> \$1160</p>	<input type="checkbox"/>	

	<p><b>Arched Foam Pad-</b> Add on Memory foam pad for knee flexion contracture. Mounts easily to back wall of knee support by Velcro. <b>EACH</b></p>	<p>A: 5" B: 5"</p>	<p><b>SB31 001</b> \$90 EACH</p>	<input type="checkbox"/>
--	---	--------------------	----------------------------------	--------------------------

	<p><b>Sandals w/Ankle Snug</b> Offers unique magnetic buckle which can be operated single-handed. <b>Pair</b></p>	<p>Size 1- A: 5" - 9" B: 5" C: 9"</p>	<p><b>S19</b> \$387</p>	<input type="checkbox"/>
	<p>Size 2- A: 7" - 12" B: 5" C: 9"</p>	<p><b>S21</b> \$438</p>	<input type="checkbox"/>	

## OPTIONAL COMPONENTS

<b>TRAY AND ACCESSORIES</b>		<p><b>Tray</b> Offers angle and depth adjustment for the perfect fit.</p> <p>A: 23" B: 16" C: 10"</p>	<p><b>STZ019</b>    \$659</p>
		<p><b>Tray Extension Wings</b> Attach to the rear of the tray to offer resting position for elbows. <b>PAIR</b> Tray Not Included</p> <p>A: 6" B: 6"</p>	<p><b>STZ035</b>    \$228</p>
<b>LOWER BODY SUPPORT</b>		<p><b>Mini Kit</b> Developed to facilitate standing for smaller children, helps to improve muscle development and weight-bearing through the feet. Gray <b>*recommend acquiring knee supports for future use when the child reaches 22"</b></p> <p>A: 7" B: 13"</p>	<p><b>STZ022</b>    \$655</p>
		<p><b>Tray Holder</b> Offers storage for the tray when not being used.</p>	<p><b>STZ010</b>    \$84</p>
		<p><b>Elbow Blocks</b> Prevents retraction of the arms on the tray support surface. <b>PAIR</b> <b>*Requires Tray Extension Wings</b></p> <p>A: 4" B: 3" C: 6" Pad Angle: 65°</p>	<p><b>UBEB</b>    \$286</p>

## OPTIONAL COMPONENTS- Continued

<b>SHOULDER SUPPORT</b>		<p><b>Flat shoulder support</b> A simple, yet effective wide pad for supporting the shoulders when using the Standz in a Supine position</p> <p>A: 4" B: 12" C: 2"</p>		<p>MH09-0** is illustrated but any of the Head Positioning aids shown can be used to control the limitation of head movement.</p>
		<p><b>STZ014-**</b>    \$254</p>		
		<p><b>Shoulder support with protractors</b> Offers angle adjustment and promotes shoulder girdle stability</p> <p>A: 5" B: 2" C: 6" D: 4" E: 2" F: 9"</p>		
	<p><b>STZ026-**</b>    \$576</p>			
<b>PELVIC AND BACK SUPPORT</b>		<p><b>Upper Back Pad</b> To control hyper extension and promote a more neutral standing position <b>**Pelvic Support Required</b></p> <p>A: 6" B: 8"</p>		<p>This configuration shows the Standz being used in Prone with the added benefit of support to the posterior of the child to prevent them from pushing into extension. Ideal for scenarios where head control is developing through prone standing but where there is a risk of the child over extending through tonal movement.</p>
	<p><b>BACK-STZ016</b></p> <p><b>PELVIC-STZ015</b></p>	<p><b>STZ016-**</b>    \$168</p>		
		<p><b>Pelvic Support Pad</b> Provides a firmer support pad and aids in achieving hip extension and neutral pelvic alignment.</p> <p>A: 6" B: 8"</p>		
	<p><b>STZ015-**</b>    \$659</p>			