



BILLING/SHIPPING INFORMATION

| Date: | Name: |
|-----------------|-----------|
| PO Number: | Address: |
| Branch: | |
| Mark for: | City: |
| Dealer Contact: | State: |
| Email Address: | Zip Code: |
| Phone: | Phone: |
| 0 | |



A versatile standing system offering supine or upright standing that accommodates abducted or neutral leg positioning for children between approx 1-9 years.





Configuration for representation only

KEY FEATURES -Tool free adjustments -Individual leg abduction and adjustment -Anti-microbial covers to aid infection control in multi-use environments -Built in gauges enable precise positioning

| Age Range (approx) | Size 1 Supine 6 months* - 9 years |
|--|---|
| User Height | 21"* - 54" |
| User Weight (max) | 99 Lbs |
| Chest/Hip Pad Width (Min/Max) | 7"- 12" |
| Height from footplates to top of upper support | 23" - 40" |
| Height of Stander (floor to top of upper pad, upright position) | 43" |
| Height of Stander (floor to support board, horizontal position) | 32" |
| Kneeblock Width (center to center, in neutral standing position) | 5" - 8" |
| Footplate Angle Range | 25 to 10 degrees |
| Individual Leg Abduction | Neutral to 30 degrees |
| Angle Range | 0 - 90 degrees |

Suggested HCPCs code E0638*

*Measurements can be achieved by the use of additional mini kit accessories (ST022)


*It is the responsibility of the CRT provider to determine the HCPCS code most appropriate for each circumstance. The HCPCS codes for reference only and are not to be deemed as legal advice. Actual coding and coverage varies by funding entities and the use of HCPCS does not guarantee coverage or payment.







CHOOSE YOUR COLOR- All part numbers containing ** will color match









REQUIRED COMPONENTS: continued

OR

| | | adjustable side wa angle and depth a Tool free adjustm | Support Width alls mounted on an adjustable bracket. hent. PAIR *Straps uired 3" - 4" |
|------------------|-------------------|--|--|
| | | oose one option below | v |
| KNEE CONTROL | | desired postion ar | trap holds knee in Id prevents flexion. AIR |
| U U | R | A: 7" B | : 17" |
| NE | | MU54-** Size 1 | |
| × | | A: 10" I | |
| | | MU68-** Size 2 | \$175 |
| | $A \rightarrow C$ | for cuffing of t control internal rotatio | nee Straps Allows he child's leg to or external knee n. PAIR |
| | | A: 15" B: 12 MU45-** Size 1 | |
| | | A: 18" B: 11 | |
| | | MU69-** Size 2 | |
| FOOT POSITIONING | | for effective leg | lid foot placement positioning and ng. Pair 3" B: 7" |
| 00 | | Size 2- A: | |
| Ĩ | - | \$2 | \$205 |

| A | Arched Foam Pad | upports- Includes and support straps i-axis bracket. PAIR |
|------|--------------------------------|---|
| | | " Max knee rence: 12" |
| COS. | STZ017-** | \$1072 |
| | | " Max knee |
| | | rence: 18" |
| | STZ023-** | \$1160 |
| | - | |
| | foam pad fo contracture. Mo | d - Add on Memory r knee flexion unts easily to back ort by Velcro. EACH |
| | A: 5" E | 3: 5" |
| AB | SB31 001 | \$89 Each |


| B | Sandals w/Ankle magnetic buck operated singl | le which ca | an be |
|---|--|----------------------------|--------------|
| | Size 1- A: 5" - 9 | ' B: 5" C: 9 | ," |
| | S19 | \$387 | |
| | Size 2- A: 7" - 12 | " B: 5" C: 9 | Э" _ |
| | S21 | \$438 | |







OPTIONAL COMPONENTS










OPTIONAL COMPONENTS





| HeadAloft includes: | |
|---------------------------------|----|
| -Adjustable, removable, | |
| padded headrest with (2) | |
| buckle attachment straps-3 | _ |
| different sized, flip-away head | |
| support straps -Occipital | 25 |
| support cushion | |

