US Retail Price List



Account Information			Satespa	ce" Hi-L
Date	Account Numb	per		
Account Name	Contact Name	9		
Phone	Fax			
Email Address				
PO Number	Mark For			
Name				e e
Address			Each Safespace Hi-l fabricated. The imag	e above is an
City	State Z	ip	example of multiple	custom selections.
Al 27 \$35,620 Safespace Hi-lo wi back wall. Viewing required below.	th standard mattress with panel in head and foot e			
ONE selection per side/end is REQU	JIRED: Main Side	Opposite Side	Head End	Foot End
Roll up wall w/ Viewing Panel (Option 1 See pg 3)	N/C	AI 28 \$1,131	N/A	N/A
Nursing Flap (Option 2 See pg 3)	AI 32 \$754	Al 28 \$1,131 Al 32 \$754	N/A	N/A
Frequent Access Opening (Option 3 See pg 3)	AI 96 \$706	AI 28 \$1,131 AI 96 \$706	N/A	N/A
Nursing Flap & Frequent Access Opening (Option 4 See pg 3)	AI 32 \$754 AI 96 \$706	Al 28 \$1,131 Al 32 \$754 Al 96 \$706	N/A	N/A
Solid wall without Viewing Panel (Option 5 See pg 3)	N/A	☐ N/C	□ N/C	N/C
Non-Opening Viewing Panel (Option 6 See pg 3)	N/A	Al 29 \$706	☐ N/C	N/C
'D' shaped Zip Opening w/ Viewing Panel	N/A	N/A	Al 31 S565	AI 31 \$565

sales@adaptive imports.com

Tel: (877) 767-9462 Fax: (877) 394-3465

(Option 7 See pg 3)









OPTIONS Cont.

Me	mory Foam Mattress	Upgrade	AI 108	1 1 4530	t suitable for individuals with ing or biting behaviors.
Fo: 30°	am Wedge and Cylino (20" tall x36" wide)(der (9" dia. x 36" wide)	AI 68	\$2,683	
	am Wedge Only (20" tall x 36" wide)		AI 69	\$1,718	
Fo: (9"	am Cylinder Only dia. x 36" wide)		AI 70	\$964	
Water Resistant Wall Seams			AI 101	\$905	
Wall		Floor	PVC Mattress	Lycra Mattress	External Corner
Light Blue		Dark Blue	Dark Blue	Black	Light Blue
Light Pink		Fuchsia	Fuchsia	*This is the top 1/2 color. Bottom 1/2 is PVC Mattress Color selected left.	Light Pink
Gray		Plum	Plum	NOT SUITABLE FOR USERS WHO HAVE BITING/ CHEWING BEHAVIORS.	Gray
Lilac		Dark Green	Dark Green		Lilac
Light Green			*If selecting Lycra Mattress Upgrade, this is the bottom 1/2 color.		Light Green
Sand					Sand
Cream		COLORS (All colors a	are approximate, variations	s to be expected)	Cream
Aqua					Aqua
					Dark Blue
					Fuchsia
					Plum




Fax: (877) 394-3465

Main Side and Opposite Side Configurations


Option 1

The zip is slightly below the top of the mattress. This is the **standard** opening for the main side and optional for the opposite side. (Al 28)




Option 2

The main zip is slightly below the top of the mattress. The nursing access panel zip is 6" above the top of the mattress. (AI 32)




Option 3

The main zip is slightly below the top of the mattress. The frequent access opening zip is 8" above the top of the mattress and opens 20" along the mattress length. (Al 96)




Option 4

The main zip is slightly below the top of the mattress. The nursing access panel zip i s6: above the top of the mattress. The frequent access opening zip is 8"" above the top fo the mattress and opens 20" along the mattress length. (Al 32 + Al 96)




Option 5

The standard configuration of the opposite side is a solid wall without window. Option 5 = solid panel for both opposite wall and Head/Foot selection.



Option 6


Non Opening viewing panel along opposite side. (AI 29)



Head and Foot End Configurations


Option 5

Optional solid wall in head or foot end for a lower stimulation environment. (No Charge)




Option 6

Standard configuration of Head and Foot end. Non opening viewing panel.



Option 7

Optional 'D' shaped zip opening for Head and/or Foot end. (Al 31)



sales@adaptiveimports.com


Tel: (877)767-9462 Fax: (877) 394-3465











Example of tube hole positioning on the main panel.

Options for tube holes: Tube holes will be made 1" diameter unless you otherwise specified. AI 30 Tubing Holes 1" each \$118 ea Which panel on are the tube Number of tube Diameter if different holes to be positioned? holes required to standard 1" Front Panel- head end Front Panel- foot end Back Panel- head end Back Panel- foot end \$118 ea AI 20 **Zip Covers** Zip covers can be sewn on one end of the zipper seam inside the compartment to prevent the occupant from accessing the backside of the zipper. Describe locations below.

sales@adaptiveimports.com Tel: (877)767-9462 Fax: (877) 394-3465



Specifications

Internal dimensions *

- Internal dimensions: 78" x 36"
- Internal height foam floor to ceiling = 50"
- Internal height top of mattress to ceiling= 44"
- Mattress thickness: 6"
- *Custom Sizes available

External dimensions

- Length x width x height: 84" x 46" x 54"
- External height at full extension Floor to top: 87"
- External height at lowest position Floor to top: 71"

Roll up Wall Opening

 Length: 62" Height: 41"

Sleeping Surface Height Adjustment

23" - 39" (Floor to Mattress Surface)

Clearance required from hard surfaces

- Total sideways travel for Hi-Lo movement: 4.5" (In head/foot direction due to lift mech. design)
- From top of external frame to room ceiling: 6"
- The front and rear opening panels need to be positioned at least 22" away from any wall or other hard surfaces
- The head and front panels need to be positioned at least 10" away from any wall or other hard surfaces
- Minimum room ceiling height required: 93"

Safe Working Load

Standard model: 220 lbs Heavy Duty: 550 lbs

> **It is the responsibility of the CRT provider to determine the HCPCS code most appropriate for each circumstance. The HCPCS codes for reference only and are not to be deemed as legal advice. Actual coding and coverage varies by funding entities and the use of HCPCS does not guarantee coverage or payment.









