

Account Information

---

Date Account Number

---

Account Name Contact Name

---

Phone Fax

---

Email Address

---

PO Number Mark For

---

Name

---

Address

---

City State Zip



















Each Safespace Siesta Hi-Lo is custom fabricated. The image above is an example of multiple custom selections.

**Safespace Siesta Hi-Lo** The Siesta is designed for those who may have a lot of movement and strength, but are unable or disinclined to pull up or climb over the high sides. The fully customizable design makes it suitable for those with complex or varying nursing needs & those requiring a mechanical lift. Its robust 'flex-out' sides can be kicked, punched, head-butted or rocked against, & help keep a person safe from injuries & from falling against hard surfaces, for example during seizures. E1399\*\*

**AI 74**  \$34,306 Safespace Siesta Hi-Lo with standard mattress with PVC cover, main zip opening w/ viewing panel. Solid back wall. Viewing panel in head and foot end, water resistant floor seams. Configuration selections required below.

<b>ONE selection per side/end is REQUIRED:</b>	<b>Main Side</b>	<b>Opposite Side</b>	<b>Head End</b>	<b>Foot End</b>
Roll up wall w/ Viewing Panel <b>(Option 1</b> See pg 3)	<input type="checkbox"/> N/C	AI 72 <input type="checkbox"/> \$1,293	N/A	N/A
Nursing Flap <b>(Option 2</b> See pg 3)	AI 32 <input type="checkbox"/> \$603	AI 72 + AI 32 \$1,293 + \$754	N/A	N/A
Frequent Access Opening <b>(Option 3</b> See pg 3)	AI 96 <input type="checkbox"/> \$565	AI 72 + AI 96 \$1,293 + \$706	N/A	N/A
Nursing Flap & Frequent Access Opening <b>(Option 4</b> See pg 3)	AI 32 + AI 96 \$603 + \$565	AI 72 + AI 32 + AI 96 \$1,293 + \$754 + \$706	N/A	N/A
Solid wall without Viewing Panel <b>(Option 5</b> See pg 3)	N/A	<input type="checkbox"/> N/C	<input type="checkbox"/> N/C	<input type="checkbox"/> N/C
Non-Opening Viewing Panel <b>(Option 6</b> See pg 3)	N/A	AI 29 <input type="checkbox"/> \$706	<input type="checkbox"/> N/C	<input type="checkbox"/> N/C
'D' shaped Zip Opening w/ Viewing Panel <b>(Option 7</b> See pg 3)	N/A	N/A	AI 73 <input type="checkbox"/> \$646	AI 73 <input type="checkbox"/> \$646

Memory Foam Mattress Upgrade	AI 108	<input type="checkbox"/>	\$539	**Not suitable for individuals with chewing or biting behaviors.	
Foam Pad covered in PVC to be attached to outside of Siesta with buckles on back for opposite side panel	AI 88	<input type="checkbox"/>	\$2,151		
Foam Pad covered in PVC to be attached to outside of Siesta with buckles on back for Head or Foot panel	AI 89	<input type="checkbox"/>	<b>Head End</b> \$1,783	<input type="checkbox"/>	<b>Foot End</b> \$1,783
Water resistant wall seams	AI 102	<input type="checkbox"/>	\$905		
Foam Wedge and Cylinder 30° (20" tall x 36" wide) (9" dia. x 36" wide)	AI 68	<input type="checkbox"/>	\$2,683		
Foam Wedge Only 30° (20" tall x 36" wide)	AI 69	<input type="checkbox"/>	\$1,718		
Foam Cylinder Only (9" dia. x 36" wide)	AI 70	<input type="checkbox"/>	\$964		

Wall	Floor	PVC Mattress	Lycra Mattress	External Corner
 <input type="checkbox"/>	 <input type="checkbox"/>	 <input type="checkbox"/>	 <input type="checkbox"/>	 <input type="checkbox"/>  <input type="checkbox"/>
 <input type="checkbox"/>	 <input type="checkbox"/>	 <input type="checkbox"/>	<p>*This is the top 1/2 color. Bottom 1/2 is PVC Mattress Color selected left. NOT SUITABLE FOR USERS WHO HAVE BITING/CHEWING BEHAVIORS.</p>	 <input type="checkbox"/>  <input type="checkbox"/>
 <input type="checkbox"/>	 <input type="checkbox"/>	 <input type="checkbox"/>		 <input type="checkbox"/>  <input type="checkbox"/>
 <input type="checkbox"/>	 <input type="checkbox"/>	 <input type="checkbox"/>		 <input type="checkbox"/>  <input type="checkbox"/>
 <input type="checkbox"/>				 <input type="checkbox"/>  <input type="checkbox"/>
 <input type="checkbox"/>				 <input type="checkbox"/>
 <input type="checkbox"/>				
 <input type="checkbox"/>				


\*If selecting Lycra Mattress Upgrade, this is the bottom 1/2 color.

**COLORS** (All colors are approximate, variations to be expected)

## Main Side and Opposite Side Configurations


### Option 1

The zip is slightly below the top of the mattress. This is the **standard** opening for the **main** side and **optional** for the **opposite side**. (AI 28)




### Option 2

The main zip is slightly below the top of the mattress. The nursing access panel zip is 6" above the top of the mattress. (AI 32)




### Option 3

The main zip is slightly below the top of the mattress. The frequent access opening zip is 8" above the top of the mattress and opens 20" along the mattress length. (AI 96)



### Option 4

The main zip is slightly below the top of the mattress. The nursing access panel zip is 6" above the top of the mattress. The frequent access opening zip is 8" above the top of the mattress and opens 20" along the mattress length. (AI 32 + AI 96)




### Option 5

The standard configuration of the opposite side is a solid wall without window. Option 5 = solid panel for both opposite wall and Head/Foot selection.



### Option 6


Non Opening viewing panel along opposite side. (AI 106 + AI 107)



## Head and Foot End Configurations


### Option 5

Optional solid wall in head or foot end for a lower stimulation environment. (No Charge)




### Option 6

Standard configuration of Head and Foot end. Non opening viewing panel.




### Option 7

Optional 'D' shaped zip opening for Head and/or Foot end. (AI 31)



**Tube Hole Options**



**Example of tube hole positioning on the main panel.**

**Options for tube holes:** Tube holes will be made 1 in. diameter unless you otherwise specified.

AI 30 Tubing Holes 1" each 1 2 3 4 5 6 \$118 ea

Which panel on are the tube holes to be positioned?	Number of tube holes required	Diameter if different to standard 1 in.
Front Panel- head end		
Front Panel- foot end		
Back Panel- head end		
Back Panel- foot end		

**Notes**

## Specifications

### Internal dimensions

- Internal dimensions: 79" x 35.5"
- Internal height from bed base to top of fabric: 55"
- Internal height from top of mattress to top of fabric: 49"
- Mattress thickness: 6"

### External dimensions

- Length x width: 88.5" x 46"
- External height at full extension: 88"
- External height at lowest position: 72.5"

### Front Opening (Long Side)

- Length: 59"
- Height: 47"

### Side Opening Door (Short Side)

- Length: 31"
- Height: 43"

### Height Adjustment

- 23" - 39" (Floor to Mattress Surface)

### Clearance required from hard surfaces

- Total sideways travel for Hi-Lo movement: 4.5" (In head/foot direction due to lift mech. design)
- From top of poles to room ceiling: 2"
- The Siesta front and rear opening panels need to be positioned at least 22" away from any wall or other hard surfaces
- The head and front panels need to be positioned at least 12" away from any wall or other hard surfaces
- Minimum floor to ceiling height required: 90"

### Safe Working Load

- Standard model: 265 lbs

\*\*It is the responsibility of the CRT provider to determine the HCPCS code most appropriate for each circumstance. The HCPCS codes for reference only and are not to be deemed as legal advice. Actual coding and coverage varies by funding entities and the use of HCPCS does not guarantee coverage or payment.