



BILLING/SHIPPING INFORMATION					
Date:		Name:			
PO Number:		Address:			
Branch:					
Mark for:		City:			
Dealer Contact:		State:			
Email Address:		Zip Code:			
Phone:		Phone:			
	Quote	Order			



Configuration for representation only

A versatile standing system offering prone or supine standing that accommodates abducted or neutral leg positioning for children between 5 - 14 years





KEY FEATURES:

-Tool free adjustments
-Highly configurable
-Individual leg abduction and adjustment
-Lift spring powered tilt function actuated by
foot pedal
Wide range of growth

	Prone	Supine
Age Range approx)	5 - 14 years	5 - 14 years
User Height	40" - 63"	40" - 63"
User Weight (max)	154 Lbs	154 Lbs
Chest Width	8.7" - 14.2"	8.7" - 14.2"
Height from footplate to top of chestplate	32.3" - 45.7"	32.3" - 45.7"
Height of Stander (floor to top of upper pad, upright position)	44.5"	48.4"
Height of Stander (floor to support board, horizontal position)	35.4"	35.4"
Kneeblock Width (center to center, in neutral standing position)	8.7"	8.7"
Footplate Angle (dorsi flexion - plantar flexion)	15° - 15°	15° - 15°
Leg Abduction	0° - 30°	0° - 30°
Angle Range	0° - 90°	0° - 90°

Suggested HCPCS code E0641**

**It is the responsibility of the CRT provider to determine the HCPCS code most appropriate for each circumstance. The HCPCS codes for reference only and are not to be deemed as legal advice. Actual coding and coverage varies by funding entities and the use of HCPCS does not guarantee coverage or payment.











CHOOSE YOUR COLOR- All part numbers containing ** will color match



Blue 01





Green 03

REQUIRED COMPONENTS



Compact Base with gas strut assisted angle adjustment and incorporated angle gauge and safety lock.

A: 51.6" B: 29.5"

\$1,631

MAST



Mast with dual leg posts and individual footplates

Pink

02

Х STZ-M2 \$1,029

OR

OR



Support Boards. In color of choice. Top board height adjustable to accommodate growth. PAIR

STZ-B2-** \$887 Х


TRUNK SUPPORT

STZ-BS2

Flip Away Laterals with Strap. Firm lateral support pads with incorporated support strap can be used to facilitate extension or rotation. PAIR

A: 3.7" B: 1.2" C: 4.7"

STZ024-** \$452




SUPPORT BOARDS

Flip Away Flexi Laterals with Strap. Flexible lateral support pads wrap around the body to provide firm yet


comfortable positioning. PAIR

A: 4.5" B: 1" C: 7.3" STZ025-** \$468



Flip Away Laterals with Strap. Firm lateral support pads with incorporated support strap can be used to facilitate extension or rotation. PAIR

A: 3.7" B: 1.2" C: 4.7" STZ024-** \$452



Flip Away Flexi Laterals with Strap.

Flexible lateral support pads wrap around the body to provide firm yet comfortable positioning. PAIR

A: 4.5" B: 1" C: 7.3"

STZ025-** \$468












REQUIRED COMPONENTS- continued














OPTIONAL COMPONENTS


















OPTIONAL COMPONENTS










Waistcoat Harness Gives full trunk and shoulder girdle stability when using the Standz 2 in the supine position

	A: 6.1" B: 5.5"					
	HB62- Size 2	\$174				
	A: 7.1" B: 6.9"					
,	HB64- Size 3	\$180				

SHOULDER SUPPORT



Shoulder support with Protractors

Shoulder protractors mounted to a flat shoulder pad, offers angle adjustment and promotes shoulder girdle stability.

A: 5.1" B: 1.6" C: 6.1" D: 3.5" E: 2" 8.7


E	F
A ((D	
C	-
В	· U

STZ026-** \$554

A: 8.5" B: 1.6" C: 5.9" D: 7.5" E: 2" F:

STZ034-** \$728

PELVIC AND BACK SUPPORT



Upper Back Support To control hyper extension and promote a more neutral standing position. *Requires Pelvic pad

A: 9.8" B: 7.7" STZ016-** \$162

Pelvic Support Pad Provides a firmer support pad and aids in achieving hip extension and neutral pelvic alignment.

A: 9.8" B: 7.7" **STZ032** \$677

This configuration shows the Standz being used in Prone with the added benefit of support to the posterior of the child to prevent them from pushing into extension. Ideal for scenarios where head control is developing through prone standing but where there is a risk of the child over extending through tonal movement.

Reset Form

Print Form

Email QUOTE Form

Email ORDER Form







F:

