US Retail Price List



| Account Information | | | Satespac Sle | :e° |
|--|--|--|--|--|
| Date | Account Numb | er | • Sie | eping |
| Account Name | Contact Name | | Comp | artmen |
| Phone | Fax | | | |
| Email Address | | | | |
| PO Number | Mark For | | | |
| Name | | | Samuel and the same of the sam | Company of the Compan |
| Address | | | | |
| City | State Zi | p | Each Safespace is cu The image above is a multiple custom select | n example of |
| Safespace Sleeping Compartment Sle | individuals who require the eping Compartment offers the face for easy transfers. Har vents any contact with hard | ne robust safety enclos d surfaces are not pado | ure of the Safespace, ded like other products | with an elevated |
| viewing panel. So | ng Compartment with stan blid back wall. Viewing par ections required below. | | | |
| CEEE Matabase to Floor Florght. | | | | |
| ONE selection per side/end is REQ | UIRED: Main Side | Opposite Side | Head End | Foot End |
| Roll up wall w/ Viewing Panel (Option 1 See pg 3) | N/C | AI 28 \$1,131 | N/A | N/A |
| Nursing Flap (Option 2 See pg 3) | Al 32 🔲 \$754 | Al 28 \$1,131 | N/A | N/A |
| | | Al 32 \$754 | | |
| Frequent Access Opening (Option 3 See pg 3) | AI 96 S565 | AI 28 \$1,131 AI 96 \$565 | N/A | N/A |
| Nursing Flap & Frequent Access Opening (Option 4 See pg 3) | AI 32 \$754 | Al 28 \$1,131 | | |
| (Option 4 See pg 5) | Al 96 \$706 | Al 32 \$754 Al 96 \$706 | N/A | N/A |
| Solid wall without Viewing Panel (Option 5 See pg 3) | N/A | N/C | N/C | N/C |
| Non-Opening Viewing Panel (Option 6 See pg 3) | N/A | AI 29 \$706 | ☐ N/C | N/C |
| 'D' shaped Zin Opening w/ Viewing Panel | N/A | N/A | Al 31 \$565 | AL 31 \$565 |

sales@adaptiveimports.com Tel: (877) 767-9462

Fax: (877) 394-3465

(Option 7 See pg 3)



| Memory Foam Mat | tress Upgrade | AI 108 | 1 1 4530 | le for individuals with iting behaviors. |
|--|--------------------------------------|--|--|--|
| Foam Wedge and 0 30° (20" tall x36" wid | Cylinder de) (9" dia. x 36" wide) | AI 68 | \$2,683 | |
| Foam Wedge Only 30° (20" tall x 36" w | ide) | AI 69 | \$1,718 | |
| Foam Cylinder Only (9" dia. x 36" wide) | y | AI 70 | \$964 | |
| Water Resistant Wa | all Seams | AI 101 | \$905 | |
| COLORS (All colors are app | proximate, variations to be expec | ted) | | |
| Wall Color (Req) | Floor Color (Req) | PVC Mattress Color (Req) | Lycra Mattress Color (Optional) | External Corner Color (Req) |
| Light Blue | Dark Blue | Dark Blue | Black | Light Blue |
| Light Pink | Fuchsia | Fuchsia | *This is the top 1/2 color. Bottom 1/2 is PVC Mattress Color selected left. | Light Pink |
| Gray | Plum | Plum | NOT SUITABLE FOR USERS WHO HAVE BITING/ CHEWING BEHAVIORS. | Gray |
| Lilac | Dark Green | Dark Green | | Lilac |
| Light Green | | *If selecting Lycra Mattress Upgrade, this is the bottom 1/2 color. | | Light Green |
| Sand | | | | Sand |
| Cream | | | | Cream |
| Aqua | | | | Aqua |
| | | | | Dark Blue |
| | | | | Fuchsia |

sales@adaptiveimports.com Tel: (877) 767-9462 Fax: (877) 394-3465




Plum

Main Side and Opposite Side Configurations


Option 1

The zip is slightly below the top of the mattress. This is the **standard** opening for the **main** side and **optional** for the **opposite side**. (Al 28)




Option 2

The main zip is slightly below the top of the mattress. The nursing access panel zip is 6" above the top of the mattress. (Al 32)




Option 3

The main zip is slightly below the top of the mattress. The frequent access opening zip is 8" above the top of the mattress and opens 20" along the mattress length. (AI 96)



Option 4

The main zip is slightly below the top of the mattress. The nursing access panel zip is 6" above the top of the mattress. The frequent access opening zip is 8" above the top of the mattress and opens 20" along the mattress length. (AI 32 + AI 96)




Option 5

The standard configuration of the opposite side is a solid wall without window.



Option 6


Non Opening viewing panel along opposite side.
(Al 29)



Head and Foot End Configurations


Option 5

Optional solid wall in head or foot end for a lower stimulation environment. (No Charge)




Option 6

Standard configuration of Head and Foot end. Non opening viewing panel.




Option 7

Optional 'D' shaped zip opening for Head and/or Foot end. (Al 31)



Tube Hole Options



Example of tube hole positioning on the main panel.

Options for tube holes: Tube holes will be made 1" diameter unless you otherwise specified.

| AI 30 | Tubing Holes 1" each 1 2 3 4 5 6 \$118 ea | | | | |
|-------|---|----------------------------------|---|--|--|
| | Which panel on are the tube holes to be positioned? | Number of tube holes required | Diameter if different to standard 1" | | |
| | Front Panel- head end | | | | |
| | Front Panel- foot end | | | | |
| | Back Panel- head end | | | | |
| | Dook Danal foot and | | | | |

| | Back Panel- foot end | | |
|-------|----------------------|---|--|
| AI 20 | Zip Covers | 1 2 3 4 5 6 | \$118 ea |
| • | | ne end of the zipper seam inside the backside of the zipper. D | de the compartment to prevent escribe locations below. |
| | | | |



Specifications

Internal dimensions *

- Internal dimensions: 78" x 36"
- Internal height foam floor to ceiling = 50"
- Internal height top of mattress to ceiling= 44"
- Mattress thickness: 6"
- *Custom Sizes available

External dimensions

18" mattress height: 84"L x 46"W x 65"H
24" mattress height: 84" L x 46"W x 71"H

Roll up Wall Opening

Length: 62"Height: 41"

Base Clearance

- 18" mattress height has 10" clear space under frame.
- 24" mattress height has 16" clear space under frame.

Clearance required from hard surfaces

- From top of external frame to room ceiling: 6"
- The front and rear opening panels need to be positioned at least 22" away from any wall or other hard surfaces
- The head and front panels need to be positioned at least 10" away from any wall or other hard surfaces
- Minimum room ceiling height required: 93"

Safe Working Load

• 396 lbs

**It is the responsibility of the CRT provider to determine the HCPCS code most appropriate for each circumstance. The HCPCS codes for reference only and are not to be deemed as legal advice. Actual coding and coverage varies by funding entities and the use of HCPCS does not guarantee coverage or payment.

