





Dreama™ by Jenx

The Dreama[™] by Jenx is a support system for effective 24hr postural management, suitable for children and adults.

Easy to use on a number of bed frames including profiling hospital beds. Allows for good ventilation and airflow

Can be used to position the user in an abducted lying position to aid hip development

- Infinitely adjustable positions thanks to the uniquely designed glide lock system
- Easy to keep clean and hygienic

Age range: 9 months to Adult Max user weight: 242lbs

Key

features

Always Supportive...



EMAIL & WEB:

sales@adaptiveimports.com www.adaptiveimports.com

PHONE & FAX: Tel: (877) 767-9462

2760 Circleport Dr. Erlanger, KY 41018 Fax: (877) 394-3465

ADDRESS:



distributors of

An Improved Quality of Life.

The Dreama[™] system by Jenx is a key part of any 24hr postural management program.

Recognizing that abnormal and posturally destructive tone is present 24/7, the Dreama[™] system follows on from support offered by conventional positioning. Utilizing a variety of available support and positioning pads, the system can be configured to support even the most complex of postures.

Supporting at risk limbs and reducing tonal spasms, the Dreama[™] contributes towards more recuperative periods of rest, which in turn has a positive contribution towards minimizing the risk of muscular skeletal deformity. Dreama[™] can be used as a preventative system or to reduce the risk of further postural deterioration, and is appropriate for any age.

Dreama is ideal for any client who has a prescribed supported seating system and looks to maximize the positive impact of this system by continuing to work against the destructive influences of abnormal tone.

*Bed frame not included

showing the hugely effective way in which Dreama prevents a destructive posture and continues all the corrective and preventative work done by conventional 'daytime' seating positions.







Custom Profiles

The flexible aluminum base contours to create the lying surface of choice for each user, whether it's achieved by placing the base of Dreama™ on a profiling hospital bed or using rolls, wedges or cushions. This means it can be used for a range of different user requirements in a variety of environments (home, care homes, support centers etc).



Individual Covers

Individual covers mean any that are soiled can be removed and replaced without the need to remove the user from their bed, minimizing sleep disruption and the need for assisted transfers.



Glide-lock[™] System

The unique Glide-lock™ system makes positioning simple! The base of the pad is inserted in to

INNOVATIVE FEATURE

Modular Cushions

INNOVATIVE FEATURE

Dreama[™] 24 hour

positioning system

Each modular support cushion is 6" deep and upholstered in fireretardant, water resistant, 2-way stretch, antibacterial fabric for infection control, as well as being breathable. The spaces between cushions allow airflow, helping to maintain a comfortable temperature. The individual cushions enable specific areas to be cleaned easily during the night with minimum disturbance.



Customizable and **Configurable Positioning**

Hundreds of individual specifications can be configured, specific to the user's needs. The extensive range of positioning options is possible because of the innovative modular system that allows the supports to slide easily into channels between the individual cushions. There are up to 16 available channels for the supports meaning that you can achieve exactly the position required.



the channel on the base and the pad is gently pushed in to place. Once in place, simply let go and pushing on the pad won't cause it to move. To remove the pad, simply hold it and carefully slide it out. There are also 'Locking Feet', which can be screwed into place within the channels (they are removable and reusable) either to mark the place a pad should go (making it easier for carers to position a user) or to put behind the pad to prevent users with strong repetitive movements from gradually 'walking' the pad out of place during the night.

Flexible Aluminum Base

The base is fixed to the bed using safety straps and the postural support cushions are then clipped into place.

Dreama[™] by Jenx Technical Specification

	Units	Cot 2006	Twin Size DC15
Age Range (approx)	mo/yrs	9 m 4 yrs	9 mo - adult
Base Size (length x width)	inch	47 22 22	78.7 x 35.4
Height Overall	inch	5.9	5.9
Single Cushion Unit size (length x width x depth)	inch	2 3. 1 5.9	35.4 x 3.7 x 5.9
No. of Cushion Units supplied	single:double	<u> </u>	16:2
Max User Weight	lbs	176	242

The Dreama™ is a modular product, allowing you to choose specific accessories to meet the exact requirements of each user.



Dreama[™] support with natural sheets Includes:

- Flexible aluminum base
- Full set of mattress cells
- Full set of neutral covers •
- Pair of anchor straps *Bed frame not included

Always Configurable Accessories

	Description Code	Dim	ensions		Neutral Covers for Dreama Accessories	
	hygienic PU and available in Large -	- DM09 4"(l) ×	3½"(w) x6½"(h)	Co Small - Mediur Large -	Il Neutral over - DU22 m - DU25 - DU28 e - DU31	Optional Memory Foam Pad Optional Cover for Memory Foam Pad Small DP01 Small DU70 Medium DP02 Medium DU73 Large DP03 Large DU76 X Large DP04 X Large DU94
R	Slim Support Pad Ideal for supporting the trunk, this slimmer version of the standard support pad allows for more comfortable arm positioning.	DM36	6¼"(h) x 4¼" (w)		Optio Cove DU96	nal Neutral Optional Padded r Cover DU95
	Abduction Block This can be used to ensure hips are abducted and the body is straight. Ideal for users with increase adductor tone to help with hip preservation.			Optional Neutral Cover Size 1 - DU52 Size 2 - DU55		
1	Knee Supports (Pair) Use these to comfortably hold knees in a specified position.	Small - DL01 Large - DL02	3"(w) x 4"(d) x 2"(h 4¼"(w) x 4"(d) x 4"		Optio Small - Large -	
÷	Height Adjustable T-Roll Offers knee support and prevents scissoring. Asymmetrically adjustable with removable pommel.	DL03	20¾"(w) x 4¼"(d) >	∢ 4¾"(h)	Optio DU61	onal Neutral Cover
	Side Lying Pad Available in different thicknesses to encourage a neutral leg separation when side lying.	Size 2 - DP06 Size 3 - DP07	17½"(w) x 5½"(d) x 21¼"(w) x 6¼ d) x		Size 2	onal Neutral Cover - DU132 - DU133
	Side Lying Support Strap Can be used to secure the upper leg when side lying. Requires use with 2 support pads.	Medium 5½"(w) x 15¼"(l) x 2¼"(h) DU141		2¼"(h)	Optional Neutral Cover Medium - DU143	
	Positioning Hoops A great functional item, to keep bedclothes away from the body or allow movement within a limited range for a comfortable night's sleep.	Small - DB01 Large - DB03	17"(l) x 2¾"(w) x 1 27½"(l) x 2¾"(w) x		Optio Small - Large -	
	Angle Adjustable Footrest Ingeniously designed to support feet in supine to prevent feet dropping.	DL05			Optio DU67	onal Neutral Cover
	Locking Feet (Set of 6) These useful items are used to mark the position where an accessory should be placed. It keeps accessories in place for users with high movement.	, DL09	5½"(w) x 15¼"(l) x	: 2¼"(h)		
	u 🔌					distributors of



EMAIL & WEB: sales@adaptiveimports.com www.adaptiveimports.com **PHONE & FAX:** Tel: (877) 767-9462 Fax: (877) 394-3465 ADDRESS: 2760 Circleport Dr. Erlanger, KY 41018 JENX Positioning for life