



Safespaces Client/Product Consideration

It is very important that you choose the right SafeSpaces product to meet the needs of the individual using it. The Adaptive Imports team can guide you through the process of choosing the appropriate product & designing it to meet your specific requirements. You can choose the size & color of your safe room or bed as well as features, such as where windows & blinds will be, or whether you want a white wall for use with projectors or other sensory equipment. The below questionnaire can be used to help direct the medical professionals towards the optimal SafeSpaces option for clients.

Patient Information and Considerations:

- 1. Medical Diagnosis-** *confirm diagnoses and medical need.*
- 2. Medical and Behavioral needs-** *sensory outbursts, self-injurious, safety, elopement, climbing/scaling, medical and ADL care.*
- 3. Level Supervision Required:** *constant 1:1, needs space for safety of self and others, elopement, poor cognition, poor balance.*
- 4. Mobility:** *fully mobile or limited mobility, full enclosure or open top, height of wall*
- 5. Sleep Patterns, Sleep Issues:** *awake at night, wonders at night.*
- 6. Intended Use of SafeSpaces Product:** *sleep, sensory space, daytime use, adapted therapy activities.*
- 7. Environment to be used:** *Home (room set up), required space measurement for SafeSpaces, needs option for travel.*

Full Enclosure Options: with Canopy top	SafeSpace (room) Sleeping Compartment Hi-Lo
Open Top Options:	Siesta Siesta Hi-Lo CosyFit CosyFit Hi-Lo CosyFit Profiling- <i>only option with adjustable bed features</i>
Travel Option-	Voyager Voyager Comact

Features and Accessories:

Standard Features	Roll up Wall with Clear PVC Viewing Panel Padded floor Foam Mattress/PVC Cover Water Resistant Floor Seams Color Options: Wall, Floor, PVC Mattress, External Corner Cover
--------------------------	--

Accessory Descriptions	**product specific accessories noted on order forms
Windows	With or without Zipper
Removeable Window Blind	Velcro PVC cover on outside to allow for darkened space
Mattress	Memory Foam/PVC cover Memory Foam/Soft-touch cover (black) Full Floor Mattress (SafeSpace room)
Postural Support	Positioning Wedge and Bolster <i>*to support postural change when profiling not option</i>
Zip Access	Access just above mattress line too allow fiber optic/sensory supports
Zip Back Covers	Inhibit opening SafeSpaces from inside
Visual Deterrent	PVC material to block space between product and wall or ceiling
Projector options	White viewing wall on back of SafeSpace(room) Ceiling PVC Porthole for projector or video camera view
Tube Holes	Ports to help control medical tubing
Floor/Wall Fixing	Bracket Hardware and Velcro option