



BILLING/SHIPPING INFORMATION

| Date: | | Name: |
|-----------------|-------|-----------|
| PO Number: | | Address: |
| Branch: | | |
| Mark for: | | City: |
| Dealer Contact: | | State: |
| Email Address: | | Zip Code: |
| Phone: | | Phone: |
| | QUOTE | ORDER |

A versatile standing system offering prone, supine or upright standing that accommodates abducted or neutral leg positioning for children between approx 1-9 years.





KEY FEATURES -Tool free adjustments -Individual leg abduction and adjustment -Anti-microbial covers to aid infection control in multi-use environments -Built in gauges enable precise positioning



| | Size 1 Prone | Size 1 Supine |
|--|-----------------------|-----------------------|
| Age Range (approx) | 6 months* - 9 years | 6 months* - 9 years |
| User Height | 21"* - 54.3" | 21"* - 54.3" |
| User Weight (max) | 99 Lbs | 99 Lbs |
| Chest/Hip Pad Width (Min/Max) | 7"- 12" | 7"- 12" |
| Height from footplates to top of upper support | 23" - 39" | 23" - 39" |
| Height of Stander (floor to top of upper pad, upright position) | 43" | 43" |
| Height of Stander (floor to support board, horizontal position) | 32" | 32" |
| Kneeblock Width (center to center, in neutral standing position) | 5" - 8" | 5" - 8" |
| Footplate Angle Range | 10 - 25 degrees | 25 to 10 degrees |
| Individual Leg Abduction | Neutral to 30 degrees | Neutral to 30 degrees |
| Angle Range | 0 - 90 degrees | 0 - 90 degrees |

Suggested HCPCs code E0641*

*Measurements can be achieved by the use of additional mini kit accessories (ST022)


*It is the responsibility of the CRT provider to determine the HCPCS code most appropriate for each circumstance. The HCPCS codes for reference only and are not to be deemed as legal advice. Actual coding and coverage varies by funding entities and the use of HCPCS does not guarantee coverage or payment.







CHOOSE YOUR COLOR- All part numbers containing ** will color match







Standz 2024 Order Form



| | | REQUIRED COMPO | DNENT | S: continued |
|------------------|-------------------|---|-------|--------------|
| | A | Standard Knee Support Width adjustable side walls mounted on an angle and depth adjustable bracket. Tool free adjustment. PAIR A: 5" B: 5" C: 3" - 4" STZ 007-** \$758 | OR | A |
| | | on below- required for Supine use | | |
| KNEE CONTROL | | Standard Knee Strap holds knee in desired postion and prevents flexion. PAIR | | |
| H | R | A: 7" B: 17" | | |
| N. | W SIR | MU54-** Size 1 \$134 | | |
| – | | A: 10" B: 21" | | A |
| | | MU68-** Size 2 \$175 | | |
| | $A \rightarrow C$ | De-Rotational Knee Straps Allows for cuffing of the child's leg to control internal or external knee rotation. PAIR | | |
| | | A: 15" B: 12" C: 12" MU45-** Size 1 \$152 | | |
| | | MU45-** Size 1 \$152 A: 18" B: 11" C: 11" | | |
| | | MU69-** Size 2 \$175 | | |
| FOOT POSITIONING | | Sandals Ensure solid foot placement for effective leg positioning and standing. Pair | | B |
| PC | | Size 1- A: 3" B: 7" | | |
| от | | S16 \$178 | | |
| 5 0 | | Size 2- A: 4" B: 9" | | |
| | | S2 \$205 | | |

| OR | A | Multi-Axis Knee Supports- Includes Arched Foam Pad <i>and</i> support straps mounted on multi-axis bracket. PAIR | | |
|----|----------------------------------|--|---------------------|--|
| | 2 | A: 5" B: 5" Max knee | | |
| | | circumference: 12" | | |
| | | STZ017-** | \$1072 | |
| | | A: 5.3" B: 4. | 5" Max knee | |
| | | circumference: 18" | | |
| | | STZ023-** | \$1160 | |
| | | | | |
| | | Arched Foam Pac | - Add on Memory | |
| | | foam pad for knee flexion | | |
| | contracture. Mounts easily to ba | | unts easily to back | |
| | | wall of knee support by Velcro. EACH A: 5" B: 5" | | |
| | | | | |
| | AB | SB31 001 | \$89 EACH | |


| B C | Sandals w/Ankle Snug Offers unique magnetic buckle which can be operated single-handed. Pair | | |
|--------|--|----------|--|
| 20 | Size 1- 5" - 9" B | 5" C: 9" | |
| | Size 2- 7" - 12" B: 5" C: 9" | | |
| Cor. | S21 | \$438 | |







OPTIONAL COMPONENTS










OPTIONAL COMPONENTS



HEADALOFT



