



BILLING/SHIPPING INFORMATION

	DIEENVO/JIIII III	
Date:		Name:
PO Number:		Address:
Branch:		
Mark for:		City:
Dealer Contact:		State:
Email Address:		Zip Code:
Phone:		Phone:
	Quote	Order

A versatile standing system offering prone or supine standing that accommodates abducted or neutral leg positioning for children between 5 - 14 years





KEY FEATURES: -Tool free adjustments -Highly configurable -Individual leg abduction and adjustment -Lift spring powered tilt function actuated by foot pedal -Wide range of growth



	Prone	Supine
Age Range approx)	5 - 14 years	5 - 14 years
User Height	40" - 63"	40" - 63"
User Weight (max)	154 Lbs	154 Lbs
Chest Width	9" - 149"	9" - 149"
Height from footplate to top of chestplate	32" - 46"	32" - 46"
Height of Stander (floor to top of upper pad, upright position)	45"	49"
Height of Stander (floor to support board, horizontal position)	35"	35"
Kneeblock Width (center to center, in neutral standing position)	9"	9"
Footplate Angle (dorsi flexion - plantar flexion)	15° - 15°	15° - 15°
Leg Abduction	0° - 30°	0° - 30°
Angle Range	0° - 90°	0° - 90°
		<i>c i i i i</i>

Suggested HCPCS code E0641**

**It is the responsibility of the CRT provider to determine the HCPCS code most appropriate for each circumstance. The HCPCS codes for reference only and are not to be deemed as legal advice. Actual coding and coverage varies by funding entities and the use of HCPCS does not guarantee coverage or payment.

sales@adaptiveimports.com Tel: (877) 767-9462 Fax: (877)394-3465



ADAP	TIVE IMPORTS (Standz 2 2024 Order Form							
	СНОО	SE YOUR COLO	R- All part nun	nbers co	ontaining '	** will color	[.] match		
	Blue 01			Pink 02				Green	
			REQUIRED C	OMPO	NENTS				
BASE	gas str angle a and inc angle		STZ-M2	post indiv	th dual leg is and vidual plates	SUPPORT BOARDS	color of c board adjust accom growt	Boards. In hoice. Top height able to modate h. PAIR	
TRUNK SUPPORT		lateral support pac support strap can	Is with Strap. Firm Is with incorporated be used to facilitate rotation. PAIR ' C: 5" \$470	OR	A B B		Flip Away Flexi La Flexible lateral so around the body to comfortable po A: 5" B: 1 STZ025-**	upport pads v to provide firr ositioning. PA	vrap m yet
PELVIC SUPPORT		lateral support pac support strap can	Is with Strap. Firm Is with incorporated be used to facilitate rotation. PAIR	OR	A B C		Flip Away Flexi La Flexible lateral so around the body to comfortable po A: 5" B: 1 STZ025-**	upport pads v to provide firr ositioning. PA	vrap m yet









0





REQUIRED COMPONENTS- continued




	R	Multi-Axis Knee Supports- Includes Arched Foam Pad <i>and</i> support straps mounted on multi-axis bracket. PAIR				
R		Size 1 (STZ017) Max strap length- 12"				
		STZ017-**	\$1072			
		Size 2 (STZ023) Max strap length- 18"				
	A: 5.3" B: 4.5"	STZ023-**	\$1160	1 [








OPTIONAL COMPONENTS











OPTIONAL COMPONENTS В 1 \sim Tray Can be used in all applications. Offers angle and depth adjustment for the perfect fit. Multi-Position- A: 27" B: 18" C: 22" STZ1028-** \$860 Prone- 27" B: 18" C: 10" TRAY STZ1029-** \$860 Tray Extension Wings Attach to the rear Elbow Blocks Prevents retraction of the of the tray to offer resting position for arms on the tray support surface. PAIR elbows. Tray Not Included *Requires Tray Extension Wings A: 6" B: 6" A: 4" B: 3" C: 6" Pad Angle: 65" UBEB STZ035 \$286 \$228 Waistcoat Harness Gives full trunk and shoulder girdle stability when using the *IHORACIC* SUPPORT Standz 2 in the supine position A: 6" B: 6" HB62- Size 2 \$181 A: 7" B: 7" HB64- Size 3 \$187 Shoulder support with Protractors Shoulder protractors mounted to a flat SHOULDER SUPPORT shoulder pad, offers angle adjustment and promotes shoulder girdle stability. This configuration A: 5" B: 2" C: 6" D: 4" E: 2" F: shows the Standz 8.7 being used in Prone with the added benefit STZ026-** \$576 of support to the A: 9" B:2" C: 6" D: 8" E: 2" F: 9 posterior of the child to prevent them from STZ034-** pushing into \$757 extension. Ideal for scenarios where head Upper Back Support: Control hyper ВАСК R control is developing STZ016 extension & promote neutral standing through prone **PELVIC AND BACK** position. *Requires Pelvic pad standing but where there is a risk of the A:10" B: 8" SUPPORT PELVIC child over extending STZ032 STZ016-** \$168 through tonal movement. Pelvic Support Pad Provides a firmer support pad and aids in achieving hip extension and neutral pelvic alignment. A: 10" B: 8" STZ032 \$704

sales@adaptiveimports.com Tel: (877) 767-9462 Fax: (877)394-3465









Additional External Head Support Option				
HeadAloft includes:				
-Adjustable, removable,				
padded headrest with (2)				
buckle attachment straps-3				
different sized, flip-away head				
support straps -Occipital	25			
support cushion				
HEADALOFT	\$490			

