

BILLING/SHIPPING INFORMATION				
Date:			Name:	
PO Number:			Address:	
Branch:				
Mark for:			City:	
Dealer Contact:			State:	
Email Address:			Zip Code:	
Phone:			Phone:	
QUOTE		<input type="checkbox"/>	ORDER	
		<input type="checkbox"/>		

A versatile standing system offering prone or upright standing that accommodates abducted or neutral leg positioning for children between 5 - 14 years



KEY FEATURES:


- Tool free adjustments
- Highly configurable
- Individual leg abduction and adjustment
- Lift spring powered tilt function actuated by foot pedal
- Wide range of growth

	Prone
Age Range approx)	5 - 14 years
User Height	40" - 63"
User Weight (max)	154 Lbs
Chest Width	9" - 14"
Height from footplate to top of chestplate	32" - 46"
Height of Stander (floor to top of upper pad, upright position)	45"
Height of Stander (floor to support board, horizontal position)	35"
Kneeblock Width (center to center, in neutral standing position)	9"
Footplate Angle (dorsi flexion - plantar flexion)	15° - 15°
Leg Abduction	0° - 30°
Angle Range	0° - 90°

Suggested HCPCS code E0638**

***It is the responsibility of the CRT provider to determine the HCPCS code most appropriate for each circumstance. The HCPCS codes for reference only and are not to be deemed as legal advice. Actual coding and coverage varies by funding entities and the use of HCPCS does not guarantee coverage or payment.*

CHOOSE YOUR COLOR- All part numbers containing ** will color match



REQUIRED COMPONENTS

BASE		Compact Base with gas strut assisted angle adjustment and incorporated angle gauge and safety lock. A: 52" B: 30"	MAST		Mast with dual leg posts and individual footplates	SUPPORT BOARDS		Support Boards. In color of choice. Top board height adjustable to accommodate growth. PAIR	
	STZ-BS2-PRN	\$1713		<input checked="" type="checkbox"/>	STZ-M2		\$1080	<input checked="" type="checkbox"/>	




TRUNK SUPPORT		Flip Away Laterals with Strap. Firm lateral support pads with incorporated support strap can be used to facilitate extension or rotation. PAIR A: 4" B: 1" C: 5"	OR		Flip Away Flexi Laterals with Strap. Flexible lateral support pads wrap around the body to provide firm yet comfortable positioning. PAIR A: 5" B: 1" C: 7"
	STZ024-**	\$470		<input type="checkbox"/>	STZ025-**

PELVIC SUPPORT		Flip Away Laterals with Strap. Firm lateral support pads with incorporated support strap can be used to facilitate extension or rotation. PAIR A: 4" B: 1" C: 5"	OR		Flip Away Flexi Laterals with Strap. Flexible lateral support pads wrap around the body to provide firm yet comfortable positioning. PAIR A: 5" B: 1" C: 7"
	STZ024-**	\$470		<input type="checkbox"/>	STZ025-**




REQUIRED COMPONENTS- continued

KNEE CONTROL		<p>Standard Knee Supports Width adjustable side walls mounted on an angle and depth adjustable bracket. Tool free adjustment. PAIR.</p> <p>A: 5" B: 5" C: 3" - 4"</p> <table border="1"> <tr> <td>STZ007-**</td> <td>\$758</td> <td><input type="checkbox"/></td> </tr> </table>	STZ007-**	\$758	<input type="checkbox"/>	OR		<p>Multi-Axis Knee Supports- Includes Arched Foam Pad and support straps mounted on multi-axis bracket. PAIR</p> <p>Size 1 (STZ017) Max strap length- 12"</p> <table border="1"> <tr> <td>STZ017-**</td> <td>\$1072</td> <td><input type="checkbox"/></td> </tr> </table> <p>Size 2 (STZ023) Max strap length- 18"</p> <table border="1"> <tr> <td>STZ023-**</td> <td>\$1160</td> <td><input type="checkbox"/></td> </tr> </table> <p>A: 5.3" B: 4.5"</p>	STZ017-**	\$1072	<input type="checkbox"/>	STZ023-**	\$1160	<input type="checkbox"/>
	STZ007-**	\$758	<input type="checkbox"/>											
	STZ017-**	\$1072	<input type="checkbox"/>											
	STZ023-**	\$1160	<input type="checkbox"/>											
		<p>Standard Knee Strap holds knee in desired position and prevents flexion. PAIR. *Optional</p> <p>A: 7" B: 17"</p> <table border="1"> <tr> <td>MU54-** Size 1</td> <td>\$134</td> <td><input type="checkbox"/></td> </tr> </table> <p>A: 10" B: 21"</p> <table border="1"> <tr> <td>MU68-** Size 2</td> <td>\$175</td> <td><input type="checkbox"/></td> </tr> </table>	MU54-** Size 1	\$134	<input type="checkbox"/>	MU68-** Size 2	\$175	<input type="checkbox"/>		<p>Arched Foam Pad- Add on Memory foam pad for knee flexion contracture. Mounts easily to back wall of knee support by Velcro. EACH</p> <p>A: 5" B: 5"</p> <table border="1"> <tr> <td>SB31 001</td> <td>\$89 EACH</td> <td><input type="checkbox"/></td> </tr> </table>	SB31 001	\$89 EACH	<input type="checkbox"/>	
	MU54-** Size 1	\$134	<input type="checkbox"/>											
MU68-** Size 2	\$175	<input type="checkbox"/>												
SB31 001	\$89 EACH	<input type="checkbox"/>												
	<p>De-Rotational Knee Straps Allows for cuffing of the child's leg to control internal or external knee rotation. PAIR *Optional</p> <p>A: 15" B: 12" C: 12"</p> <table border="1"> <tr> <td>MU45-** Size 1</td> <td>\$152</td> <td><input type="checkbox"/></td> </tr> </table> <p>A: 18" B: 11" C: 11"</p> <table border="1"> <tr> <td>MU69-** Size 2</td> <td>\$175</td> <td><input type="checkbox"/></td> </tr> </table>	MU45-** Size 1	\$152	<input type="checkbox"/>	MU69-** Size 2	\$175	<input type="checkbox"/>							
MU45-** Size 1	\$152	<input type="checkbox"/>												
MU69-** Size 2	\$175	<input type="checkbox"/>												
FOOT POSITIONING		<p>Sandals Ensure solid foot placement for effective leg positioning and standing. Pair</p> <p>Size 2- A: 4" B: 9"</p> <table border="1"> <tr> <td>S2</td> <td>\$203</td> <td><input type="checkbox"/></td> </tr> </table> <p>Size 3- A: 4.1" B: 10.8"</p> <table border="1"> <tr> <td>S3</td> <td>\$218</td> <td><input type="checkbox"/></td> </tr> </table>	S2	\$203	<input type="checkbox"/>	S3	\$218	<input type="checkbox"/>		<p>Sandals w/Ankle Snug Offers unique magnetic buckle which can be operated single-handed. Pair</p> <p>Size 2- :7" - 12" B: 5" C: 9"</p> <table border="1"> <tr> <td>S21</td> <td>\$438</td> <td><input type="checkbox"/></td> </tr> </table> <p>Size 3- N/A</p>	S21	\$438	<input type="checkbox"/>	
	S2	\$203	<input type="checkbox"/>											
	S3	\$218	<input type="checkbox"/>											
S21	\$438	<input type="checkbox"/>												

OPTIONAL COMPONENTS

TRAY		<p>Tray Can be used in all applications. Offers angle and depth adjustment for the perfect fit.</p> <p>A: 27" B: 18" C: 10"</p> <p>STZ1029-** \$860</p>	
		<p>Tray Extension Wings Attach to the rear of the tray to offer resting position for elbows. Tray Not Included</p> <p>A: 6" B: 6"</p> <p>STZ035 \$228</p>	

OPTIONAL COMPONENTS

SHOULDER SUPPORT		<p>Shoulder support with Protractors Shoulder protractors mounted to a flat shoulder pad, offers angle adjustment and promotes shoulder girdle stability.</p> <p>A: 5" B: 2" C: 6" D: 4" E: 2" F: 9"</p> <p>STZ026-** \$576</p>		<p><i>This configuration shows the Standz being used in Prone with the added benefit of support to the posterior of the child to prevent them from pushing into extension. Ideal for scenarios where head control is developing through prone standing but where there is a risk of the child over extending through tonal movement.</i></p>
		<p>A: 9" B: 2" C: 6" D: 8" E: 2" F: 9"</p> <p>STZ034-** \$757</p>		
PELVIC AND BACK SUPPORT	 <p>BACK STZ016</p>	<p>Upper Back Support To control hyper extension and promote a more neutral standing position.</p> <p><i>*Requires Pelvic pad</i></p> <p>A: 10" B: 8"</p> <p>STZ016-** \$168</p>		
	 <p>PELVIC STZ032</p>	<p>Pelvic Support Pad Provides a firmer support pad and aids in achieving hip extension and neutral pelvic alignment.</p> <p>A: 10" B: 8"</p> <p>STZ032 \$704</p>		