

BILLING/SHIPPING INFORMATION					
Date:			Name:		
PO Number:			Address:		
Branch:					
Mark for:			City:		
Dealer Contact:			State:		
Email Address:			Zip Code:		
Phone:			Phone:		
		<input type="checkbox"/>	QUOTE	<input type="checkbox"/>	
		<input type="checkbox"/>	ORDER	<input type="checkbox"/>	

A versatile standing system offering supine or upright standing that accommodates abducted or neutral leg positioning for children between approx 1-9 years.



Configuration for representation only

KEY FEATURES

- Tool free adjustments
- Individual leg abduction and adjustment
- Anti-microbial covers to aid infection control in multi-use environments
- Built in gauges enable precise positioning

	Size 1 Supine
Age Range (approx)	6 months* - 9 years
User Height	21"* - 54.3"
User Weight (max)	99 Lbs
Chest/Hip Pad Width (Min/Max)	6.9" - 12"
Height from footplates to top of upper support	22.8" - 39.4"
Height of Stander (floor to top of upper pad, upright position)	42.5"
Height of Stander (floor to support board, horizontal position)	31.9"
Kneeblock Width (center to center, in neutral standing position)	5.1" - 8.3"
Footplate Angle Range	25 to 10 degrees
Individual Leg Abduction	Neutral to 30 degrees
Angle Range	0 - 90 degrees

*Measurements can be achieved by the use of additional accessories (ST022)

Suggested HCPCS code E0638*

**It is the responsibility of the CRT provider to determine the HCPCS code most appropriate for each circumstance. The HCPCS codes for reference only and are not to be deemed as legal advice. Actual coding and coverage varies by funding entities and the use of HCPCS does not guarantee coverage or payment.*

CHOOSE YOUR COLOR- All part numbers containing ** will color match

Blue

01		
-----------	--	--

Pink

02		
-----------	--	--

Green

03		
-----------	--	--

Red

04		
-----------	--	--

REQUIRED COMPONENTS

BASE		Compact base with gas strut assisted angle adjustment and incorporated angle gauge and safety lock. A: 35" B: 25.4"	MAST		Mast with dual leg posts and individual footplates	SUPPORT BOARDS		Support Boards. Pair in color of choice. Top board height adjustable to accommodate growth.
	STZ-BS1-SUP	\$1,495		<input checked="" type="checkbox"/>	STZ-M1		\$947	<input checked="" type="checkbox"/>

TRUNK SUPPORT		Lateral Supports with Rotation Control Strap Width and height adjustable - can be set to accommodate asymmetry. With incorporated rotation control strap for simple and effective control and stability. A: 3.7" B: 1.2" C: 4.7"
		STZ012-** \$448

OR

	Rotation Control Strap A wide support belt without lateral pads, threaded through the support fixings. Fastened securely around child and then pulled to facilitate extension or rotation. A: 5" B: 19.7"
	STZ013-** \$448

PELVIC SUPPORT		Lateral Supports with Rotation Control Strap Width adjustable with incorporated rotation control strap for simple and effective control and stability. A: 3.7" B: 1.2" C: 4.7"
		STZ012-** \$448

OR

	Rotation Control Strap A wide support belt without lateral pads, threaded through the support fixings. Fastened securely around child and then pulled to facilitate extension or rotation. A: 5" B: 19.7"
	STZ013-** \$448

REQUIRED COMPONENTS: continued

KNEE CONTROL		<p>Standard Knee Support Width adjustable side walls mounted on an angle and depth adjustable bracket. Tool free adjustment. PAIR *Straps required</p>	<p>A: 4.5" B: 4.9" C: 2.6" - 3.7"</p> <p>STZ 007-** \$729</p>	OR		<p>Multi-Axis Knee Supports- Includes a domed knee pad and support straps mounted on multi-axis bracket. PAIR</p>	<p>A: 5.3" B: 4.5" Max knee circumference: 12.2"</p> <p>STZ017-** \$1,031</p>
	*Please choose one option below				<p>A: 5.3" B: 4.5" Max knee circumference: 17.7"</p> <p>STZ023-** \$1,115</p>		
		<p>Standard Knee Strap holds knee in desired position and prevents flexion. PAIR</p>	<p>A: 6.7" B: 16.9"</p> <p>MU54-** Size 1 \$128</p> <p>A: 9.6" B: 21.1"</p> <p>MU68-** Size 2 \$167</p>			<p>Arched Foam Pad- Memory foam pad for knee flexion contracture. Mounts easily to back wall of knee support by Velcro. EACH</p>	<p>A: 4.5" B: 5.4"</p> <p>SB31 001 \$86</p>
		<p>De-Rotational Knee Straps Allows for cuffing of the child's leg to control internal or external knee rotation. PAIR</p>	<p>A: 15.4" B: 11.8" C: 11.8"</p> <p>MU45-** Size 1 \$145</p> <p>A: 18.1" B: 11.4" C: 11.4"</p> <p>MU69-** Size 2 \$167</p>				
FOOT POSITIONING		<p>Sandals Ensure solid foot placement for effective leg positioning and standing. Pair</p>	<p>Size 1- A: 3.1" B: 6.7"</p> <p>S16 \$171</p> <p>Size 2- A: 3.7" B: 8.7"</p> <p>S2 \$195</p>		<p>Sandals w/Ankle Snug Offers unique magnetic buckle which can be operated single-handed. Pair</p>	<p>Size 1- 4.7" - 8.9" B: 4.7" C: 8.9"</p> <p>S19 \$372</p> <p>Size 2- 6.5" - 11.6" B: 4.7" C: 8.9"</p> <p>S21 \$421</p>	

OPTIONAL COMPONENTS

HEAD SUPPORT		Molded Head Support		Flat Head Support																
HEAD SUPPORT	<p>Multigrip Headrest Provides height, depth and angle adjustment. Includes removable black fabric cover. Available with 1 (A+B) 2 (A+B+C) or 3 (A+B+C+D) link sets. Black</p> <p>A: 5.1" B: 3.1" C: 3.1" D: 3.1"</p> <table border="1"> <tr> <td>MH06 (A-B)</td> <td>\$681</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td>MH01 01 (A-C)</td> <td>\$740</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	MH06 (A-B)	\$681	<input type="checkbox"/>	<input type="checkbox"/>	MH01 01 (A-C)	\$740	<input type="checkbox"/>	<input type="checkbox"/>	<p>Anatomically molded PU headrest for minimal – moderate head support. Bracket for height depth and angle adjustment</p> <p>A: 8.3" B: 11"</p> <table border="1"> <tr> <td>MH07-**</td> <td>\$565</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	MH07-**	\$565	<input type="checkbox"/>	<input type="checkbox"/>	<p>For children who require minimal head positioning</p> <p>A: 3.9" B: 7.5"</p> <table border="1"> <tr> <td>MH09-**</td> <td>\$185</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	MH09-**	\$185	<input type="checkbox"/>	<input type="checkbox"/>	
	MH06 (A-B)	\$681	<input type="checkbox"/>	<input type="checkbox"/>																
	MH01 01 (A-C)	\$740	<input type="checkbox"/>	<input type="checkbox"/>																
	MH07-**	\$565	<input type="checkbox"/>	<input type="checkbox"/>																
MH09-**	\$185	<input type="checkbox"/>	<input type="checkbox"/>																	
TRAY AND ACCESSORIES	<p>Tray Can be used in all applications. Offers angle and depth adjustment for the perfect fit.</p> <p>A: 23" B: 15.6" C: 17.7"</p> <table border="1"> <tr> <td>STZ008-**</td> <td>\$634</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	STZ008-**	\$634	<input type="checkbox"/>	<input type="checkbox"/>	<p>Tray Holder Offers storage for the tray when not being used.</p> <table border="1"> <tr> <td>STZ010</td> <td>\$81</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	STZ010	\$81	<input type="checkbox"/>	<input type="checkbox"/>	<p>Tray Extension Wings Attach to the rear of the tray to offer resting position for elbows. PAIR Tray Not Included</p> <p>A: 5.9" B: 5.9"</p> <table border="1"> <tr> <td>STZ035</td> <td>\$219</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	STZ035	\$219	<input type="checkbox"/>	<input type="checkbox"/>	<p>Prevents retraction of the arms on the tray support surface.</p> <p>PAIR</p> <p>*Requires Tray Extension Wings</p> <p>A: 3.9" B: 3.3" C: 5.5" Pad Angle: 65°</p> <table border="1"> <tr> <td>UBEB</td> <td>\$275</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	UBEB	\$275	<input type="checkbox"/>	<input type="checkbox"/>
	STZ008-**	\$634	<input type="checkbox"/>	<input type="checkbox"/>																
	STZ010	\$81	<input type="checkbox"/>	<input type="checkbox"/>																
	STZ035	\$219	<input type="checkbox"/>	<input type="checkbox"/>																
UBEB	\$275	<input type="checkbox"/>	<input type="checkbox"/>																	
LOWER BODY SUPPORT	<p>Developed to facilitate standing for smaller children, helps to improve muscle development and weight-bearing through the feet.</p> <p>Gray *we recommend acquiring knee supports for future use when the child reaches 22"</p> <p>A: 6.7" B: 12.6"</p> <table border="1"> <tr> <td>STZ022</td> <td>\$630</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	STZ022	\$630	<input type="checkbox"/>	<input type="checkbox"/>															
	STZ022	\$630	<input type="checkbox"/>	<input type="checkbox"/>																

OPTIONAL COMPONENTS

SHOULDER SUPPORT		<p>Flat shoulder support A simple, yet effective wide pad for supporting the shoulders when using the Standz in a Supine position</p>
		<p>A: 3.7" B: 12.4" C: 2.4"</p>
	<p>STZ014-** \$244 <input type="checkbox"/></p>	
		<p>Shoulder support with protractors Offers angle adjustment and promotes shoulder girdle stability</p>
<p>A: 5.1" B: 1.6" C: 6.1" D: 3.5" E: 2" F: 8.7"</p>		
<p>STZ026-** \$554 <input type="checkbox"/></p>		

Reset Form

Print Form

Email QUOTE Form

Email ORDER Form